

5mm Disposable Clip Applier



Counter window

Features:

- + Big in-jaw rotation force to help prevent clip dislodgement during procedure.
- + Unique Clip and jaw design to provide secure fixation on the vessel and resistance to migration.
- + Clear counter window shows the quantity of clip remaining inside the shaft.
- + The safety mechanism prevents the empty jaws closing to damage the tissue.
- + 360° rotation shaft to facilitate visualization and accurate placement.
- + Ergonomic handle piece design for convenient manipulation.

Code No.	Item	Spec.(mm)	Units/ Box
CIT-ML0515	5mm Disposable clip applier	5 x 330	1