

5mm Disposable Clip Applier **New**



Features:

- + Smooth and reliable clip loading and firing mechanism.
- + Fired clip with good retention force on target, reduce risk of clip dislodgment.
- + 360° rotation via rotation knob, facilitate visualization and accurate placement.
- + Matted stainless steel shaft, eliminate the reflected glare from scope.
- + Empty jaw with safety mechanism, prevent undesired closure.
- + Ergonomic handle design, provide convenient manipulation.

Code No.	Item	Spec.(mm)	Units/ Box
CIT-ML0515	5mm Disposable clip applier	5 x 330	1