

AW Channel Cleaning Adapter

For Fujifilm Endoscopes 700 Series



Features

- Compatible with Fujifilm® Endoscopes 700 series.
- Bedside cleaning directly after the procedure, preventing the formation of biofilm.
- Single use only, eliminate manual cleaning and reprocessing of the reusable adapters.
- Particular Orange & Black color appearance and warning tag to eliminate misuse during endoscopic procedure.

Ref. Code	Description	Units/Box
CIT-AW-002	AW Channel Cleaning Adapter, Non-sterile	100