

# AW Channel Cleaning Adapter and Suction Cleaning Adapter



AW channel cleaning adapter



Suction cleaning adapter

## Intended Use

- AW Channel Cleaning Adapter is attached to the air/water cylinder of the endoscope, used for the air water channel pre-cleaning after the procedure.
- Suction Cleaning Adapter is used to aspirate reprocessing fluids through the instrument channel port of the endoscope.

## Features

- Bedside cleaning directly after the procedure, preventing the formation of biofilm.
- Single use only, eliminate manual cleaning and reprocessing of the reusable adapters.
- Adapters come packaged in a set for ease of use.
- Particular Yellow & Black color appearance and warning tag to eliminate misuse during endoscopic procedure.
- Compatible with Olympus® Endoscope 140/160/180/190/240/260/290 Series.

Ref. Code	Description	Units/Box
CIT-AW-001	AW Channel Cleaning Adapter, Non-sterile	100
CIT-BC-001	Suction Cleaning Adapter, 0.5 Meter Length Tube, Non-sterile	25
CIT-BC-001-EX	Extended Suction Cleaning Adapter, 1.0 Meter Length Tube, Non-sterile	25
CIT-AK-1	Kit: AW Channel Cleaning Adapter + Suction Cleaning Adapter, Non-sterile	25
CIT-AK-1-EX	Kit: AW Channel Cleaning Adapter + Extended Suction Cleaning Adapter, Non-sterile	25

