


Bariatric Surgery Stent

Application: Placed after a gastric sleeve surgery, prevents the possible staple line leak, promotes the healing of fistula, establishes the digestive tract to allow the introduction of food.

Performance Characteristics:

- Fully covered with PU/silicone
- OTW delivery system for distal release
- Extraction threads on both ends for easy removal and repositioning
- Larger diameter and long length of the stent: prevent migration
- Golden radiopaque markers for clearer visibility: 4+2+4, total 10pcs

Stent Shape Y




2. Proximal cup end, distal ball end

4. Two cup ends

8. Two ball ends



Specifications:

| Product Code | Dia.A (mm) | Total Length B (mm) | Stent Shape Y | Covering | Introducer |
|--------------|---------------|------------------------|------------------|---|----------------|
| NTJB- ABY-FC | 22/24 | 200/220/240 | 2/4/8 |  | 24Fr*700mm OTW |

Remark: FC means full covered.