


# Biliary Stent (PTCD/ERCP)


**Application:** Biliary stent is used for bile duct stenosis caused by various malignant lesions.




**Performance Characteristics:**

- Self expanding with high radial force
- Excellent positional stability
- Fully covered with PU/silicone, with good resistant and elastic performance
- Atraumatic ends
- Extraction threads on both ends for easy removal and repositioning
- 6 Radiopaque gold markers for better visibility

**Stent Shape Y**



**Specifications:**

| Type | Product Code  | Dia.A<br>(mm) | Total Length B<br>(mm)                    | Stent Shape<br>Y | Product<br>Illustration                                                               | Introducer   |
|------|---------------|---------------|-------------------------------------------|------------------|---------------------------------------------------------------------------------------|--------------|
| PTCD | NTJA3-ABY     | 8/10          | 40/50/60/70/80/90/100/110/120<br>60/70/80 | 1<br>15          |  | 8Fr*500mm    |
| ERCP | NTJA4-ABY     | 8/10          | 40/50/60/70/80/90/100/110/120<br>60/70/80 | 1<br>15          |  | 8/9Fr*1850mm |
| PTCD | NTJA3-ABY- HC | 8/10          | 40/50/60/70/80/90/100/110/120<br>60/70/80 | 1<br>15          |  | 8Fr*500mm    |
| ERCP | NTJA4-ABY- HC | 8/10          | 40/50/60/70/80/90/100/110/120<br>60/70/80 | 1<br>15          |  | 8/9Fr*1850mm |
| PTCD | NTJA3-ABY- FC | 8/10          | 40/50/60/70/80/90/100/110/120<br>60/70/80 | 1<br>15          |  | 8Fr*500mm    |
| ERCP | NTJA4-ABY- FC | 8/10          | 40/50/60/70/80/90/100/110/120<br>60/70/80 | 1<br>15          |  | 8/9Fr*1850mm |

Remark: HC means half covered; FC means full covered