

Bouffant Cap & Beard Cover



Bouffant Cap



Beard Cover

Intended use:

Ideal for use in hospitals and the care area for nurses and patients. Wherever hair and beard coverage is required, including sectors such as food and catering.

Features:

- + Made of soft and lightweight polypropylene non-woven material.
- + With double stitched rubber spring.
- + Highest wearing comfort with good fit.
- + Breathable for the best possible heat regulation.
- + Hypo-allergenic.

Code No.	Product Name	Size	Color	Raw Material
CIT-TM-21	Bouffant Cap	21"	Blue/White	PP Non-woven fabric 10±2gsm
CIT-TM-24	Bouffant Cap	24"	Blue/White	PP Non-woven fabric 10±2gsm
CIT-TM-28	Bouffant Cap	28"	Blue/White	PP Non-woven fabric 10±2gsm
CIT-HXT-18	Beard Cover	18"	Blue/White	PP Non-woven fabric 10±2gsm

Packaging:

Product	Pcs/Bag	Bags/CTN	Pcs/CTN	Carton Size(cm)
Bouffant Cap	100pcs	20bags	2000pcs	65x25x24
Beard Cover	100pcs	20bags	2000pcs	61x25x25

• Products can be customized. Please contact us for any further information.

