

Bougie Dilator Set



Intended use

- Dilation treatment of digestive tract and cardia.

Features

- Made of medical macromolecule, with proper rigidity.
- Equipped with radiopaque mark at the tip, easy to observe under X-ray.
- Equipped with nitinol super-elastic guide wire, bite block, flusher syringe.

Product code	Balloon Size (mm)	Length(mm)	Characteristics	Units/Box
CIT-BDS-6	5, 7, 9, 11, 13, 15	850	Each set is equipped with 6 or 8 bougies, 1 guide wire,	1
CIT-BDS-8	5, 7, 9, 11, 13, 15,17,19	850	1 bite block and 1 flusher syringe.	1