

Clean Cap

For Olympus



For Pentax



For Fujifilm



Features

- Deliver Air or CO2 for insufflation and water for lens cleaning.
- Compatible with Olympus®, Fujifilm® 500/600 Series, and Pentax® 90 Series Endoscopes.
- Universal bottle cap compatible for most water bottles on the market, like B.Braun, Fresenius, Baxter.
- 24 hours use, discard daily.

Product Code	Description	Units/Box
CIT-IT-001	Clean cap for Olympus Endoscopes 140/160/180/190/240/260/290 Series. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-001-WO	Clean cap without two-way valve for Olympus Endoscopes 140/160/180/190/240/260/290 Series. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-001-EX	Extended Clean cap for Olympus Endoscopes 140/160/180/190/240/260/290 Series: length of the two-way valves tube connecting is 1.1 meters. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-002	Clean cap for Fujifilm Endoscopes 500/600 series. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-002-WO	Clean cap without two-way valve for Fujifilm Endoscopes 500/600 series. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-002-EX	Extended Clean cap for Fujifilm Endoscopes 500/600 series: length of the two-way valves tube connecting is 1.1 meters. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-003	Clean cap for Pentax 90 Series Endoscopes. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-003-WO	Clean cap without two-way valve for Pentax 90 Series Endoscopes. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-003-EX	Extended Clean cap for Pentax 90 Series Endoscopes: length of the two-way valves tube connecting is 1.1 meters. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-004	Clean cap with direct connector for Pentax 90 Series Endoscopes. Compatible with B.Braun and Fresenius sterile water bottles.	10
CIT-IT-009	Clean Cap for Olympus Endoscopes 140/160/180/190/240/260/290 Series. Compatible with Asia-Pacific Baxter sterile water bottles	10