

Clean Cap for Fujifilm Endoscopes 700 Series



Features

- Deliver air or CO2 for insufflation and water for lens cleaning.
- Compatible with Fujifilm® 700 Series Endoscopes.
- Universal bottle cap compatible for most water bottles on the market, like B.Braun, Fresenius, Baxter.
- Bottle Cap for Asian& Pacific Baxter water bottle is also available.
- 24 hours use, discard daily.

Product Code	Description	Units/Box
CIT-IT-FJ700	Clean Cap for Fujifilm Endoscopes 700 Series, Compatible with B.Braun and Fresenius water bottles.	10
CIT-IT-FJ700-EX	Extended Clean Cap for Fujifilm Endoscopes 700 Series :length of the two-way valves tube connecting is 1.1 meters, Compatible with B.Braun and Fresenius water bottles.	10
CIT-IT-FJ700-KT	Clean Cap with CO2 Source Tube (for Fujifilm Endoscopes 700 Series), Compatible with B.Braun and Fresenius water bottles	15
CIT-IT-FJ700B	Clean Cap for Fujifilm Endoscope 700 Series, Compatible with Asian& Pacific Baxter water bottles.	10
CIT-IT-FJ700B-EX	Extended Clean Cap for Fujifilm Endoscope 700 Series :length of the two-way valves tube 1.1 meters, Compatible with Asian& Pacific Baxter water bottles.	10
CIT-IT-FJ700B-KT	Clean Cap with CO2 Source Tube (for Fujifilm Endoscope 700 Series), Compatible with Asian& Pacific Baxter water bottles	15