

Clean Cap with CO2 Source Tube



Metal Head with Slot



CO₂ Filter

Product Code	Description	Units/Box
CIT-IT-001-KT	Clean Cap with CO ₂ Source Tube (for Olympus Endoscopes 140/160/180/190/240/260/290 Series). Compatible with B.Braun and Fresenius sterile water bottles.	15
CIT-IT-002-KT	Clean Cap with CO ₂ Source Tube (for Fujifilm Endoscopes 500/600 Series). Compatible with B.Braun and Fresenius sterile water bottles.	15
CIT-IT-003-KT	Clean Cap with CO ₂ Source Tube (for Pentax 90 Series Endoscopes). Compatible with B.Braun and Fresenius sterile water bottles.	15
CIT-IT-009-KT	Clean Cap with CO ₂ Source Tube (for Olympus Endoscopes 140/160/180/190/240/260/290 Series). Compatible with Asia-Pacific Baxter sterile water bottles.	15

CO₂ Source Tube

Intended use

- To be used with a CO₂ insufflator to supply CO₂ to a GI endoscope during GI endoscopic procedures.

Specification

- With 0.2µm filtering membrane. Tube length: 2.15 meters.



Product Code	Description	Units/Box
CIT-CFT-001	CO ₂ Source Tube	15

