










# Cleaning Brush



### Intended use

- Cleaning of endoscope working and valve channels.

Picture	Ref. Code	Diameter of Brush (mm)	Length of Brush (mm)	Characteristics	Units/Box
	CIT-CB-V1	10	35	Endoscope Valve	100
	CIT-CB-V2	5/10	20/30	Endoscope Valve	100
	CIT-CB-0202-1800-C	2/2	20/20	≥2.0mm Endoscope Channel	100
	CIT-CB-0303-2300-C	3/3	20/20	≥2.0mm Endoscope Channel	100
	CIT-CB-0505-2300-C	5/5	20/20	≥2.8mm Endoscope Channel	100
	CIT-CB-0510-2300--C	5/10	20/30	≥2.8mm Endoscope Channel	100
	CIT-CB--KIT 1	10 + 5/5	35 + 20/20	Endoscope Valve + ≥2.8mm Endoscope Channel	100
	CIT-CB-KIT 2	5/10 + 5/5	20/30 + 20/20	Endoscope Valve + ≥2.8mm Endoscope Channel	100