


# Colon/Rectum Anastomosis Stent

**Application:** To cover colonic and rectal fistula or in case of spontaneous or iatrogenic ruptures.

**Performance Characteristics:**

- Self expanding with high radial force
- Excellent positional stability
- Half covered with PU/silicone, with good resistant and elastic performance
- Atraumatic ends
- Extraction threads on both ends for easy removal and repositioning
- 6 Radiopaque gold markers for better visibility



**Specifications:**

Product Code	Dia.A (mm)	Total Length B (mm)	Product Illustration	Introducer
NTJA8-AB-HC	20	60		18Fr*1400mm/ 24Fr*700mm OTW
	30			18Fr*650mm Proximal Release

Remark: HC means half covered.