



Disposable Auto Linear Stapler



1. Anvil
2. Protector
3. Retaining pin
4. Cartridge
5. Cartridge holder
6. Anvil holder
7. Active handle
8. Fixed handle
9. Retaining pin thumb button
10. Release button



Features:

- + Thickness of tissue sutured adjustable during closing the instrument.
- + Bigger opening distance, good for tissue arrangement.
- + Clamp tissue automatically, easy for one-hand manipulation.
- + Different kinds of staple height for thin/thick tissue.

Code No.	Cartridge article No	Suture thickness(mm)	Cartridge Color	Legs Height(mm)
ZFB30	FCB30	1.5	Blue	3.8
ZFH30	FCH30	2.0	Green	4.5
ZFB45	FCB45	1.5	Blue	3.8
ZFH45	FCH45	2.0	Green	4.5
ZFB60	FCB60	1.5	Blue	3.8
ZFH60	FCH60	2.0	Green	4.5