

Disposable Simple Retrieval Bag



Features:

- + High strength TPU material of the bag ensures no leakage.
- + Transparent bag provides a better visibility of the specimen.
- + The string helps to close and retrieve the bag easily. Short and long two kinds of string are available for choosing.
- + Two versions, U shape bag or V shape bag available for choosing.

Code No.	Size of Bag(mm)	Bag Volume (ml)	Units/Box
HKQWD-U1	100 x 135 (U Shape)	200	100
HKQWD-U2	125 x 130 (U Shape)	250	100
HKQWD-U3	160 x 165 (U Shape)	400	100
HKQWD-V1	125 x 185 (V Shape)	200	100
HKQWD-V2	165 x 210 (V Shape)	320	100
HKQWD-V3	195 x 200 (V Shape)	420	100