


## Disposable TPU Retrieval Bag



### Features:

- + High strength TPU material of the bag ensures no leakage.
- + Self opening pouch design provides a safe and easy capture of the specimen.
- + Transparent bag provides a better visibility of the specimen.
- + Simple delivery to aid rapid entry of the product.
- + Two models are available to meet various demands.

Product Model	Code No.	Size of bag(mm)	Bag volume (ml)	Diameter of Introducer(mm)	Units/ Box
Disposable TPU Retrieval Bag	CIT-RB-01	120 x 150	300	10	10
	CIT-RB-02	150 x 200	500	10	10
	CIT-RB-03	180 x 200	1000	12	10
Disposable TPU Retrieval Bag (String Type)	CIT-RB-01-S	120 x 150	300	10	10
	CIT-RB-02-S	150 x 200	500	10	10
	CIT-RB-03-S	180 x 200	1000	12	10