

Disposable Balloon Trocars



Features:

- + Universal seal, compatible for a wider size range of instruments.
- + Locking mechanism between obturator and cannula, provides a safe insertion.
- + Balloon and fixation design provides great stabilization. No damage to the abdominal wall.
- + Three kinds of tip available: Bladeless, Bladed and Optical.

Code No.	Product Name	Specification (mm)	Units/ Box
CIT-BT-05	Disposable Balloon Trocar, Optical	5 x 100	8
CIT-BT-10		10 x 100	6
CIT-BT-12		12 x 100	6
CIT-BTW-05	Disposable Balloon Trocar, Bladeless	5 x 100	8
CIT-BTW-10		10 x 100	6
CIT-BTW-12		12 x 100	6
CIT-BTY-05	Disposable Balloon Trocar, Shield Bladed	5 x 100	8
CIT-BTY-10		10 x 100	6
CIT-BTY-12		12 x 100	6