




Disposable Biopsy Forceps



Features

- Unique edge sharpness processing on the oval cup.
- Alligator jaws to prevent slippage.
- Clean tissue sampling and minimize the bleeding.
- Three length markers on the 2300mm sheath offer precise instrument location even used under gastroscopy.
- Color-coded sheath and slider for easy recognition of area compatibility.
- Two different handle designs, both with rotatable thumb ring for easy manipulation.

Ref. Code No.	Cup Type	Specifications	With Spike	Color of Sheath and Slider	Endoscope Channel	Area Compatibility
DH-BF-181050C	Oval	1.8*1050 mm	No	Blue 	2.0 mm	Bronchoscopy
DH-BF-181200C	Oval	1.8*1200 mm	No	Blue 	2.0 mm	Bronchoscopy
DH-BF-181600C	Oval	1.8*1600 mm	No	Blue 	2.0 mm	Bronchoscopy
DH-BF-231600C	Oval	2.3*1600 mm	No	Yellow 	2.8 mm	Gastroscopy
DH-BF-231600CS	Oval	2.3*1600 mm	Yes	Yellow 	2.8 mm	Gastroscopy
DH-BF-232300C	Oval	2.3*2300 mm	No	Purple 	2.8 mm	Enteroscopy
DH-BF-232300CS	Oval	2.3*2300 mm	Yes	Purple 	2.8 mm	Enteroscopy
DH-BE-231600C	Alligator	2.3*1600 mm	No	Yellow 	2.8 mm	Gastroscopy
DH-BE-231600CS	Alligator	2.3*1600 mm	Yes	Yellow 	2.8 mm	Gastroscopy
DH-BE-232300C	Alligator	2.3*2300 mm	No	Purple 	2.8 mm	Enteroscopy
DH-BE-232300CS	Alligator	2.3*2300 mm	Yes	Purple 	2.8 mm	Enteroscopy