

Disposable Bipolar Forceps with Cutting



Features:

- + Four functions in one instrument: Grasping, Dissecting, Cutting and Coagulation.
- + Bipolar power blade makes the cutting clean and precise every time.
- + Unique jaws opening design, good for vessel closure up to 4mm.
- + Bipolar electrosurgical cut with coagulation along surface of cutting blade ensures less bleeding and much safer.
- + Ergonomic handle design for comfortable grip and manipulation.
- + Compatible with all standard bipolar electrosurgical units.

Code No.	Item	Spec. (mm)	Units/ Box
CIT-BH-03	Bipolar Forceps with Cutting	5 x 330	1