

Disposable Bipolar Instruments



Integrated plug



Split plug



Features:

- + Electrosurgical current restricted to just the tissue between the jaws. The risk of patient burns is reduced significantly.
- + Better control over the area being targeted and help prevent damage to other sensitive tissues.
- + Unique cable attachment design, no interruption for the procedure.
- + Rotation knob design ensure the jaws to rotate 360° freely.
- + Ergonomic handle design to provide comfort gripping and manipulation.

Code No.	Product	Spec.(mm)	Units/Box
CIT-BH-01	Disposable Bipolar Curved Dissector, double Action	5 x 330	12
CIT-BH-02	Disposable Bipolar Fenestrated Atraumatic Grasper, double Action	5 x 330	12