

# Disposable Bite Block



CIT-KD-C1



CIT-KD-B1



CIT-KD-B2

## Features

- Reinforced dental retention rim, strong and reliable.
- Oxygenating bite block delivers oxygen to both the nose and mouth during endoscopy.
- Two kinds of straps for choice: knitted elastic strap or rubber strap.
- Both kinds of straps are latex free, low allergy.
- Size: 66Fr.
- Huge improvement on patient comfort.

Product Code	Description	Units/Box
CIT-KD-C1	With oxygen inlet, knitted elastic strap, suitable for adult	50
CIT-KD-B1	With knitted elastic strap, suitable for adult	100
CIT-KD-B2	With rubber strap, suitable for adult	100