


Disposable Circular Stapler



Code No.	Code No.(Tilt Top Anvil)	Colour of Protector	Quantity of Staples	O.D.(mm)	Cutting Dia.(mm)	Staple Height(mm)
WH-Y-17	/	White	16	17	10	3.8
WH-Y-21	/	Yellow	20	21	14	4.5
WH-Y-24	WH-Y-T-24	Grey	18	24	15	4.5
WH-Y-26	WH-Y-T-26	Red	20	26	17	4.5
WH-Y-29	WH-Y-T-29	Blue	24	29	20	4.8
WH-Y-32	WH-Y-T-32	Green	26	32	22	5.0
WH-Y-34	WH-Y-T-34	Purple	28	34	24	5.0

Features:

- + Low profile anvil,easy to insert into bowel lumen.
- + Non-slippery wing nut provides increased conform and control.
- + Audible and tactile sound feedback as an indication that the firing sequence is complete.
- + Tissue gap between anvil and cartridge,avoid tissue damage.
- + 7 sizes available,better choice for surgeons.