

10mm Disposable Clip Applier



Features:

- + Jaws angled 20° to facilitate visualization of anatomic features.
- + 20 pcs of Medium/large size titanium ligation clips are pre-loaded in the device.
- + The transparent shaft provides line of sight to the remaining clips.
- + The safety mechanism prevents the handle from being squeezed when last clip is fired and the empty jaws will not close to damage the tissue.
- + The shaft can rotate 360° smoothly by rotating the rotation knob, to facilitate visualization and accurate placement.
- + Ergonomic handle piece design for convenient manipulation.

Code No.	Item	Spec.(mm)	Units/ Box
CIT-ML1020	10mm Disposable clip applier	10 x 290	1