


## Disposable Cold Snares



### Intended Use

This product is intended to resect diminutive polyp(1-5mm) or small polyp(6-9mm) in digestive tract during endoscopy procedure.

### Features

- Cold snare technique helps to eliminate the thermal injury caused by hot snare.
- Reduce the possibility of delayed hemorrhage which is caused by hot snare's polyp resection.
- Loop made of seven braided thin 0.24mm stainless steel wires, provide clean and precise cut.
- One to one rotation between the handle and loop, provide precise positioning and increase the efficiency of procedure.

Product Code	O.D. of Outer Tube (mm) D	Working Length (mm) L	Loop Width (mm) W	Shape	Units/Box
CIT-CSD-DLW-A	1.8, 2.3	1200, 1600, 1800, 2300	10,15,20	Oval	10
CIT-CSD-DLW-C	1.8, 2.3	1200, 1600, 1800, 2300	10,15,20	Hexagonal	10