


Disposable Double-handle Auto Linear Stapler



1. Anvil
2. Cartridge
3. Retaining pin
4. Firing handle
5. Closure handle
6. Fixed handle
7. Release button



Code No.	Cartridge Model	Thickness(mm)	Cartridge Color	Staple Height(mm)
SZH-30-2	SCW-30-2.5	1.0	White	2.5
SZH-30-3	SCB-30-3.8	1.5	Blue	3.8
SZH-30-4	SCH-30-4.5	2.0	Green	4.5
SZH-60-3	SCB-60-3.8	1.5	Blue	3.8
SZH-60-4	SCH-60-4.5	2.0	Green	4.5

Features:

- + Parallel closure, uniform pressure of tissue.
- + Retaining pin avoid tissue over-flow.
- + Clamp tissue automatically, easy for one-hand manipulation.
- + Three kinds of staple height for vascular, thin and thick tissue.
- + 3-row stapleline white cartridge 30mm for vascular use available.