



With String



Features:

- + High strength TPU material of the bag ensures no leakage.
- + Transparent bag provides a better visibility of the specimen.
- + Preloaded memory wire automatically opens the bag upon deployment.
- + The string helps to close and retrieve the bag easily (EB-200, EB-400).
- + After deployment, the introducer can be withdrawn, free the cannula.
- + New U shape bag design, disperse the air pressure inside the bottom of bag when pulling out the bag from the body.

Code No.	Introducer Diameter(mm)	Size of Bag(mm)	Bag Volume (ml)	Units/Box
BBD-TSG-01	10	95 x 195	200	10
BBD-TSG-02	10	130 x 195	400	10
BBD-TSG-03	10	160 x 200	800	10
EB-200 (with String)	10	95 x 195	200	10
EB-400 (with String)	10	130 x 195	400	10