


# Disposable Endoscopic Circular Stapler



Code No.	Color of protect	Number of staples	Outer Dia. (mm)	Incision Dia. (mm)	Staple Height(mm)
ZQH17	White	16	17	10	3.8
ZQH21	Yellow	20	21	14	4.5
ZQH24	Grey	18	24	15	4.5
ZQH26	Red	20	26	17	4.5
ZQH29	Blue	24	29	20	4.8
ZQH32	Green	26	32	22	5.0
ZQH34	Purple	28	34	24	5.0

## Features:

- + Low profile anvil,easy to insert into bowel lumen.
- + Good air seal for endoscopic procedure.
- + Audible and tactile sound feedback as an indication that the firing sequence is complete.
- + Rubber inlay and Non-slippery wing nut provides increased conform and control.
- + Wisdom window indicator,no need for manual adjustment.