


Disposable Endoscopic Hemoclip



Features

- Maximum clip span to 16mm, suitable for large size tissue.
- Different working lengths available, meet different clinical needs.
- Repeatedly opening and closing of the clip via the control of the handle.
- 360° synchronously rotation performance, offering precise placement of the clip.

Product Code	Working Length (mm)	Diameter of Outer Tub (mm)	Max. Grasp Degree (°)	Max. Opening Size of Clip W (mm)	Endoscope Channel (mm)	Units/Box
CIT-ZXJ-1250-9	1250	2.3	135°	9	≥2.8	10
CIT-ZXJ-1650-9	1650	2.3	135°	9	≥2.8	10
CIT-ZXJ-1950-9	1950	2.3	135°	9	≥2.8	10
CIT-ZXJ-2350-9	2350	2.3	135°	9	≥2.8	10
CIT-ZXJ-1250-11	1250	2.3	135°	11	≥2.8	10
CIT-ZXJ-1650-11	1650	2.3	135°	11	≥2.8	10
CIT-ZXJ-1950-11	1950	2.3	135°	11	≥2.8	10
CIT-ZXJ-2350-11	2350	2.3	135°	11	≥2.8	10
CIT-ZXJ-1250-13	1250	2.5	135°	13	≥2.8	10
CIT-ZXJ-1650-13	1650	2.3	135°	13	≥2.8	10
CIT-ZXJ-1950-13	1950	2.3	135°	13	≥2.8	10
CIT-ZXJ-2350-13	2350	2.3	135°	13	≥2.8	10
CIT-ZXJ-1250-16	1250	2.3	135°	16	≥2.8	10
CIT-ZXJ-1650-16	1650	2.3	135°	16	≥2.8	10
CIT-ZXJ-1950-16	1950	2.3	135°	16	≥2.8	10
CIT-ZXJ-2350-16	2350	2.3	135°	16	≥2.8	10