


Disposable Grasping Forceps




Forceps with rat teeth

- Suit for grasping foreign body in complex shape (e.g. retrieval of non-regular shaped stent)



Forceps with prongs

- Suit for grasping excised polyp, disease tissues and etc.



Extended alligator and rat tooth Forceps

- Suit for grasping flat foreign body (e.g. coin)



Product Code	Item	Specifications	O.D.(mm)	Units/Box
CIT-FG-23-160-A1	Alligator Tooth Forceps	Working Length:160cm	2.3	10
CIT-FG-23-230-A1	Alligator Tooth Forceps	Working Length:230cm	2.3	10
CIT-FG-23-160-A2	Rat Tooth Forceps	Working Length:160cm	2.3	10
CIT-FG-23-230-A2	Rat Tooth Forceps	Working Length:230cm	2.3	10
CIT-FG-23-160-A3	Extended alligator and rat tooth Forceps	Working Length:160cm	2.3	10
CIT-FG-23-230-A3	Extended alligator and rat tooth Forceps	Working Length:230cm	2.3	10
CIT-FG-24-160-B1	Foreign Body Holding Forceps	3 prongs, Working Length:160cm	2.4	10
CIT-FG-24-230-B1	Foreign Body Holding Forceps	3 prongs,Working Length:230cm	2.4	10
CIT-FG-24-160-B2	Foreign Body Holding Forceps	4 prongs, Working Length:160cm	2.4	10
CIT-FG-24-230-B2	Foreign Body Holding Forceps	4 prongs,Working Length:230cm	2.4	10