

## Disposable Grasping Forceps (With Net)



### Features

- Smooth net deployment and retraction from the catheter.
- One to one rotation, offering precise position.
- Strong and durable nets in big size range.
- Ideal for esophageal food impactions, foreign body retrieval, tissue collection and polyp removal.

Product Code	O.D. (mm)	Working Length (mm)	Net Size (W*L) mm	Endoscope Channel (mm)	Units/Box
CIT-YWQ-TY-18x1600x25	1.8	1600	15*25	≥2.0	10
CIT-YWQ-TY-23x1600x25	2.3	1600	25*45	≥2.8	10
CIT-YWQ-TY-23x1600x35	2.3	1600	35*60	≥2.8	10
CIT-YWQ-TY-23x1600x45	2.3	1600	45*80	≥2.8	10
CIT-YWQ-TY-23x2300x25	2.3	2300	25*45	≥2.8	10
CIT-YWQ-TY-23x2300x35	2.3	2300	35*60	≥2.8	10
CIT-YWQ-TY-23x2300x45	2.3	2300	45*80	≥2.8	10