


# Disposable Hemorrhoids Stapler & STARR Stapler



1. Anvil    2. Suture hole    3. Safety    4. Wing nut



P02 Anorectal introducer



P03 Anorectal dilator



P04 Purse-string suture anoscope



P05 Suture threader

Code No.	Specification(mm)	Number of Staples(piece)	Staple Height(mm)
GH-Y-32	Φ 32	30	3.8
GH-Y-34	Φ 34	32	3.8
GH-Y-36	Φ 36	34	3.8

## STARR Stapler

Code No.	Specification(mm)	Number of Staples(piece)	Staple Height(mm)
GH-ST-32	Φ 32	30	4.5
GH-ST-34	Φ 34	32	4.5
GH-ST-36	Φ 36	34	4.5

## Features:

- + Rational arrangement of staples to reduce bleeding after procedure.
- + Tissue gap between anvil and cartridge, avoid tissue damage.
- + Big volume of cartridge housing to have enough space for prolapsed tissue.
- + Four suture holes design for easy operation.
- + Non-slippery wing nut provides increased conform and control.
- + Transparent purse-string suture anoscope with scale for better vision.
- + Two kind of dilators for surgeons' better choice.
- + Detachable or fixed anvil for surgeon's better choice.
- + Three sizes available, more options for surgeons.