



Disposable Linear Cutter



1. Cartridge
2. Anvil
3. Protector
4. Cartridge holder
5. Fixed handle
6. Active handle
7. Firing button

Code No.	Cartridge Model	Cutting Length(mm)	Suturing Length(mm)	Suturing Thickness(mm)	Staple Height(mm)	Cartridge Color
QHB55-3.8	QCB55-3.8	55	57	1.5	3.8	Blue
QHH55-4.5	QCH55-4.5	55	57	2.0	4.5	Green
QHB75-3.8	QCB75-3.8	75	77	1.5	3.8	Blue
QHH75-4.5	QCH75-4.5	75	77	2.0	4.5	Green
QHB100-3.8	QCB100-3.8	100	102	1.5	3.8	Blue
QHH100-4.5	QCH100-4.5	100	102	2.0	4.5	Green

Features:

- + One instrument with blue/green cartridges.
- + Reloadable instrument to change up to 8 cartridges for single patient.
- + Safety mechanism avoids the firing with empty cartridge
- + Full range of sizes for different procedure requirements.