

# Disposable Mini Monopolar Instruments



## Intended use:

Disposable Mini Monopolar Instruments are designed to cut, grasp, dissect and coagulate selected tissue in laparoscopic surgery. The instrument is connected through a suitable monopolar cable with the monopolar output of the electro-surgical generator.

## Features:

- + Compatible with 3mm trocar port, make the trauma to patient as small as possible.
- + Ergonomic handles offer comfort grip.
- + Light weight reduces the doctors' load during the long time surgery.
- + Colorful knobs for easier distinguish.

Code No.	Item	Spec. (mm)	Color of rotation knob	Units/ Box
CIT-XMH-01	Disposable Curved Scissors	3 x 280	Orange	12
CIT-XMH-03	Disposable Curved Dissector	3 x 280	Blue	12
CIT-XMH-04	Disposable EndoClinch with Ratchet	3 x 280	Dark Grey	12
CIT-XMH-05	Disposable Fenestrated Atraumatic Grasper with Ratchet	3 x 280	Dark Grey	12
CIT-XMH-08	Disposable Blunt with Ratchet	3 x 280	Dark Grey	12
CIT-XMH-09	Disposable Fenestrated Long Grasper with Ratchet	3 x 280	Dark Grey	12

## Different Jaws:



Curved Scissor



Curved Dissector



Endoclinch



Fenestrated Grasper



Blunt Grasper



Fenestrated Long Grasper