

## Disposable Nylon Retrieval Bag



### Features:

- + Ultra-durable Ripstop Nylon bag, ensures no leakage and highly tear resistance performance.
- + Preloaded memory wire automatically opens bag upon deployment.
- + Flexible metal ring provide enough holding force to keep the bag mouth open and stable for optimal specimen loading.
- + The pulling string enables the full closure of the bag.
- + Extended introducer makes the distance not an issue anymore.

Code No.	Size of Bag (mm)	Diameter of Introducer (mm)	Working Length (mm)	Bag Volume (ml)	Units/Box
QD-A-1	80 x 100	10	290	150	12
QD-A-2	90 x 190	10	290	300	12
QD-A-3	100 x 100	10	290	175	12
QD-A-4	110 x 150	10	290	240	12
QD-A-5	110 x 200	10	290	435	12
QD-A-6	170 x 240	12	320	1100	12
QD-A-7	200 x 260	15	400	1400	12