


# Disposable Polypectomy Snare



## Features

- Made of high-strength medical braided stainless steel wire, offering a precise & quick cutting properties, as well as an outstanding electric coagulation.
- Loop rotates synchronously by rotating the 3-ring handle, highly increase efficiency.
- Standard power connection, compatible with all main high-frequency devices on the market.
- Ergonomic design of 3-ring handle, easier to hold and use.

## Oval snare

Product code	O.D.(mm)(D)	Length(mm)(L)	Loop Width(mm)(W)	Characteristics	Units/Box
CIT-SD-DLW-A1	1.8	1600	20	Standard Rotation	10
CIT-SD-DLW-A1	2.3	1600, 2300	10,15,20,25,30,35,50		10
CIT-SD-DLW-A1R	2.3	2300	20,30	Synchronous Rotation	10

## Crescent Snare

Product code	O.D.(mm)(D)	Length(mm)(L)	Loop Width(mm)(W)	Characteristics	Units/Box
CIT-SD-DLW-B1	1.8	1600	20	Standard Rotation	10
CIT-SD-DLW-B1	2.3	1600,2300	10,15,20,25,30		10

## Hexagonal Snare

Product code	O.D.(mm)(D)	Length(mm)(L)	Loop Width(mm)(W)	Characteristics	Units/Box
CIT-SD-DLW-C1	1.8	1600	20	Standard Rotation	10
CIT-SD-DLW-C1	2.3	1600,2300	10,15,20,25,30		10
CIT-SD-DLW-C1R	2.3	2300	20,30	Synchronous Rotation	10

## Cold Snare

Product code	O.D.(mm)	Length(mm)	Loop Width(mm)	Endoscope Channel(m m)	Units/Box
CIT-SD-23-230-10-C3	2.3	2300	10	≥2.8	10
CIT-SD-23-230-15-C3	2.3	2300	15	≥2.8	10

## Features

- Special designed loop with thinner wire, enables a clean cut
- Suitable for diminutive polypectomy