


Disposable Endoscopic Retrieval Bag



Disposable Endoscopic Retrieval Bag



Intended use:

Disposable Endoscopic Retrieval Bag is indicated for using in endoscopic procedures to allow tissues or debris to be removed from the abdominal cavity.

Rec. Code No.	Working Length of Introducer (mm)	Diameter of Introducer (mm)	Bag Volume (ml)	Units/ Box
QD-A-0240	290	10	240	12
QD-A-0750	320	12	750	12
QD-A-1200	320	12	1200	12
QD-A-2000	400	15	2000	12

Rec. Code No.	Recommended Anatomy Removal
QD-A-0240	Appendix, Gallbladder (Standard and Acute), Small Bowel Resection
QD-A-0750/ QD-A-1200	Bowel, Kidney, Spleen Resection, Gastric Sleeve
QD-A-2000	Enlarged Spleen, Malignant Kidney Resection, Hemi-Colectomy