

# Disposable Stone Extraction Balloon



## Intended use

- Stone removal in biliary tract, including biliary sludge and residual gallstone left in biliary tract after Mechanical Lithotripsy.

## Features

- Innovational design of three lumens (lumens for guide wire, injection and inflation) makes the insertion and radiography easier and more convenient.
- Three different sizes for each balloon, meets demands of different patients.
- Equipped with Radiopaque marks on both ends of balloon, convenient for locating under X-ray.
- Unique design of tapering catheter, improves the insertion into papilla.

Product code	O.D.(Fr)	Length(mm)	Endoscope Channel (mm)	Balloon Size(mm)	Units/Box
CIT-BAL-23-8.5/11.5/15-B	7	2000	≥ 3.2	8.5/11.5/15	1
CIT-BAL-23-9/12/15-B	7	2000	≥ 3.2	9/12/15	1
CIT-BAL-23-10/13/16-B	7	2000	≥ 3.2	10/13/16	1