



Disposable Stone Extraction Basket



Intended use

- Removal of gallstone in biliary duct and foreign bodies in upper and lower digestive tract

Features

- Innovative handle design, with functions of push, pull and rotation, easier to grasp gallstone and foreign body.
- Convenient for injection of contrast medium with injection port on handle.
- Made of advanced alloyed material, ensures good shape retention even after difficult stone removal.

Product code	O.D(mm)	Length(mm)	Basket Width (mm)	Endoscope Channel	Characteristics	Units/Box
CIT-BAS-18-70-15-N4-B-O	1.8	700	15	≥2.0	Nitinol, 4 wires, spiral shape, single-hand operation	10
CIT-BAS-18-70-20-N4-B-O	1.8	700	20	≥2.0		10
CIT-BAS-23-200-20-N4-B-D	2.3	2000	20	≥2.8	Nitinol, 4 wires, spiral shape, both hands operation	10
CIT-BAS-23-200-25-N4-B-D	2.3	2000	25	≥2.8		10
CIT-BAS-23-200-30-N4-B-D	2.3	2000	30	≥2.8	Nitinol, 4 wires, semi-circle shape, single-hand operation	10
CIT-BAS-18-70-15-N4-C-O	1.8	700	15	≥2.0		10
CIT-BAS-18-70-20-N4-C-O	1.8	700	20	≥2.0	Nitinol, 4 wires, semi-circle shape, both hands operation	10
CIT-BAS-23-200-20-N4-C-D	2.3	2000	20	≥2.8		10
CIT-BAS-23-200-25-N4-C-D	2.3	2000	25	≥2.8	Nitinol, 4 wires, semi-circle shape, both hands operation	10
CIT-BAS-23-200-30-N4-C-D	2.3	2000	30	≥2.8		10