

Disposable Surgical Face Mask



Greater safety in the operating theater and the ward

Reinforced nose clip provides seamless fit

Comfortable wear

Loops with great elasticity, less pressure on the ears

Low breathing resistance

Bigger facial covering



Disposable Surgical Face Mask

- + EN ISO13485: 2016.
- + 93/42/EEC medical product class I.
- + EN 14683: 2019 + AC: 2019, Type II / IIR.
- + Bacterial filtration efficiency(BFE) \geq 98 %.
- + 3-layer.



Composition

CITEC disposable surgical mask is composed of 3 layers: outer non-woven fabric, middle melt blown fabric and inner non-woven fabric material.

Usage

For standard and long-term use and all operating theater applications, ward use for the protection of both users and patients. A high level of comfort is guaranteed, even after prolonged use.

Ref. Code	Folded Size	Color	Package	Carton	Carton Size
CIT-FMII-L	17.5cm*9.5cm	Blue	50pcs	2000pcs	40*55*37cm



 **Changzhou Intl. Trade & Enterprises Cooperative Co., Ltd.**
 Add: No.20 Floor, Huajing Building, No.1 Hejing Garden,Zhonglou District, 213000,
 Changzhou Jiangsu, China
 Website: www.citec-group.com
 E-mail: citec@citec-group.com
 Tel: +86 519 86803321

 **Riomavix S.L.**
 Add.: Calle de Almansa 55, 1D, Madrid 28039 Spain
 Tel.: +34 658 396 230