

Disposable Trocars for Pediatric Laparoscopic Surgery



Features:

- + Low profile and less weight, specially designed for pediatric laparoscopic surgery.
- + Stainless steel cannula with thin wall and good rigidity, greatly reduce the incision and the damage to the abdominal wall.
- + Unique silicone fixator design, provide great retention force when instrument goes in and out.
- + Flexible and durable seal maintain abdominal insufflation while accommodating instrument's smooth insertion and extraction.

Code No.	Specification(mm)	Units/ Box
CIT-PT-03	3 x 50	10
CIT-PT-05	5 x 80	10