

# Disposable Wound Protector



## Features:

- + Easy and rapid set-up.
- + Provide maximum exposure with a minimum incision size.
- + Provide clear visualization and wide working space.
- + Greatly reduce the risk of wound infections.
- + Maintain moisture at the incision site.
- + Evenly dispersed tension force greatly reduce the potential risk of trauma and pain.
- + Multiple sizes available meet for different clinical demands.

Ref Code No.	Outer Ring Diameter (mm)	Inner Ring Diameter (mm)	Unfolded Channel Height (mm)	Suitable Incision Size (mm)	Units/Box
CIT-QKB60/70	60	70	150	20-40	10
CIT-QKB80/90	80	90	150	25-60	10
CIT-QKB120/130	120	130	250	50-90	10
CIT-QKB150/160	150	160	250	70-120	10
CIT-QKB180/190	180	190	250	90-140	10
CIT-QKB220/230	220	230	250	110-170	10
CIT-QKB270/280	270	280	250	170-220	10
CIT-QKB320/330	320	330	250	220-280	10