



Duodenal/Colonic/Rectal Stent with OTW Introducer System

Application: Duodenal/Rectal/Colonic stent is used for for the intestinal (duodenal,rectal,colonic) stenosis, obstruction or anastomotic stenosis.






Performance Characteristics:

- Atraumatic ends ensure patient feel less irritation in the sensitive area. Greatly improving the quality of life
- Golden radiopaque markers for clearer visibility
- Delivery system with markers for better position accuracy



Specifications:

Type	Product Code	Dia.A (mm)	Total Length B (mm)	Product Illustration	Introducer
Duodenal	NTJA7-AB-HC	18/20/22	60/80/100/120/140/160		18Fr*1400mm
	NTJA7-AB-FC				
Colonic/ Rectal	NTJA7-AB-HC	25/28/30	60/80/100/120		24Fr*700mm

Remark: HC means half covered; FC means full covered

CE 0197