


Duodenal/Colonic Stent with TTS Introducer System


Application: Duodenal/Colonic stent is used for for the intestinal (duodenal, colonic) stenosis, obstruction or anastomotic stenosis.



Performance Characteristics:

- Self-expandable and high radial force
- Reinforced metal tube for fluent release
- Ball ends for excellent positional stability
- Release under endoscopic observation with high radiopacity
- Golden radiopaque markers for clearer visibility

Specifications:

Product Code	Dia.A (mm)	Total Length B (mm)	Product Illustration	Introducer
NTJA12-AB	18/20/22/25	60/80/90/100/120/140/150		10Fr*2300mm
NTJA12-AB-HC				

Remark: HC means half covered.