


ERCP Biliary Stent

Application: The stent is used for bile duct stenosis caused by various malignant lesions as well as for the closure of postoperative or tumor-related biliary fistulas.

Performance Characteristics:

- Nitinol mesh ensures excellent expansion and flexibility of the stent.
- Silicone and PU covering available for selection.
- High radial force maintains duct patency while adapting to curved anatomy.
- Gold radiopaque markers for clearer visibility.
- Atraumatic ends minimize tissue trauma.
- Extraction lasso enables safe and easy removal of the stent.
- Capability of re-position after 80% deployment of the stent.
- Guidewire compatible up to 0.035".



Stent					Introducer System	
Product Code	A: Diameter (mm)	B: Length (mm)	Shape	Covering	Catheter Diameter(mm)	Working Length(mm)
NTJA4-ABY-HC/FC	6/8/10	30/40/50/60/70/80/90/100/120	1: Vertical Cylinder, 15: Two Flare	Uncovered/ Half Covered/ Full Covered	8	1850

Remark: HC means half covered, FC means full covered.