

# Disposable Endoscope Valves System

## For Fujifilm Endoscopes 700 Series



Air/Water valve



Suction valve






Biopsy valve



Water jet connector

### Features

- Compatible with Fujifilm® Endoscopes 700 series.
- Single use, eliminate the work burden of reprocessing.
- No lubricant on the valves to interfere the procedure or obstruct the view.
- Suction valve with metal shaft, to help maintain the suction force like the reusable valves.
- Valves come packaged in sterile kits for ease use.
- High quality and cost effective.

Product Picture	Ref. Code	Description	Units/Box
	CIT-VK-003	4-piece Valve kit (Air/Water Valve, Suction Valve, Water Jet Connector, Biopsy Valve)	25
	CIT-VK-004	3-piece Valve kit (Air/Water Valve, Suction Valve, Biopsy Valve)	25
	CIT-VK-008	2-piece Valve kit (Air/Water Valve, Suction Valve)	25

