

Disposable Endoscope Valves System For Olympus Endoscopes



Air/Water valve



Suction valve



Biopsy valve



Water jet connector (Metal Insert)










Water jet connector (All Plastic)



AW channel cleaning adapter

Features

- Compatible with Olympus® Endoscopes 140/160/180/190/240/260/290 series.
- Single use, eliminate the work burden of reprocessing.
- No lubricant on the valves to interfere the procedure or obstruct the view.
- Suction valve with metal shaft, to help maintain the suction force like the reusable valves.
- Water jet connector with unique threading design, to minimize potential to cross thread and avoid leakage.
- Valves come packaged in sterile kits for ease use.
- Great blister tray design enables the valves to drop out of tray easily.
- High quality and cost effective.

Product Picture	Ref. Code	Description	Units/Box
	CIT-VK-001	4-piece Valve kit (Air/Water Valve, Suction Valve, Water Jet Connector (Metal Insert), Biopsy Valve)	25
	CIT-VK-009	4-piece Valve kit (Air/Water Valve, Suction Valve, Water Jet Connector (All Plastic), Biopsy Valve)	25
	CIT-VK-002	3-piece Valve kit (Air/Water Valve, Suction Valve, Biopsy Valve)	25
	CIT-VK-005	2-piece Valve kit (Air/Water Valve, Suction Valve)	25
	CIT-VK-006	4-piece Valve kit (Air/Water Valve, Suction Valve, Biopsy Valve, AW Channel Cleaning Adapter)	25
	CIT-VK-007	5-piece Valve kit (Air/Water Valve, Suction Valve, Water Jet Connector (Metal Insert), Biopsy Valve, AW Channel Cleaning Adapter)	25
	CIT-VK-010	5-piece Valve kit (Air/Water Valve, Suction Valve, Water Jet Connector (All Plastic), Biopsy Valve, AW Channel Cleaning Adapter)	25