



**Quality First
Develop Harmoniously**

Rev. No.:2025A

EN ISO 13485:2016

CE

FDA



CITEC™  endopack™


ENDOSCOPY PROCEDURE KIT

Date of issue: 2025-01-17

endopack™ Endoscopy Procedure Kit

PROTECT
COLLECT
CLEAN



CITEC™  endopack™ Endoscopy Procedure Kit provides essential supplies needed for the procedure set up, bedside pre-cleaning, transport and reprocessing of flexible endoscope. It helps to improve the procedural efficiency and patient care.

Composition



Sponge with Concentrated Enzyme

Ergonomically round shape, highly match with the insertion part of endoscope. With multi enzyme detergents, to remove harmful bio-burden on the endoscopes. Sponge without concentrated enzyme is also available.



Eco Bowl

Ladder design, made with degradable material. Maximum 700ml for detergent solution.



Transport Pad

Strong leak-proof liner to prevent leakage of fluids. Cinching cords converts pad into a transport container.



Biohazard Label

Stick the Biohazard Label on the transport pad, to warn the used endoscope.




Gauze Pads

Contain 10pcs, Size: 10*10cm, 4-ply.



Lubricant Jelly

20g or 42g, two kinds for selection.



Brush Kit

Cleaning brush for Channel Opening: double ends brush, Ø5mmx40mm and Ø12mmx40mm. working length 190mm. Cleaning brush for Endoscope Channel: double tapered ends brush, Ø5mm~Ø7mmx20mm and Ø5mm~Ø7mmx20mm, tube diameter 1.7mm, working length 2300mm.



Endoscope Valves System

For Olympus endoscopes series 140/160/180/190/240/260/290. 5-piece Valve kit (Air/Water Valve, Suction Valve, Water Jet Connector, Biopsy Valve, AW Channel Cleaning Adapter).



Pre-cleaning Valve Kit

For Olympus endoscopes series 140/160/180/190/240/260/290. AW Channel Cleaning Adapter*1, Suction Cleaning Adapter*1.



Endoscope Valves System

For Olympus endoscopes series 140/160/180/190/240/260/290. 4-piece Valve kit (Air/Water Valve, Suction Valve, Water Jet Connector, Biopsy Valve).

Ref. Code	Sponge with Concentrated Enzyme	Sponge without Concentrated Enzyme	Eco Bowl 700ml	Transport Pad	Biohazard Label	Gauze Pads 10pcs	Brush Kit for Endoscope Cleaning	Lubricant Jelly, 42g	Lubricant Jelly, 20g	5-piece Valve Kit	Kits /Box
CIT-SP-001	•										50
CIT-SP-002		•									50
CIT-EK-02	•		•								10
CIT-EK-02A		•	•								10
CIT-EK-04	•		•	•	•						10
CIT-EK-04A		•	•	•	•						10
CIT-EK-06	•		•	•	•	•	•				10
CIT-EK-06A		•	•	•	•	•	•				10
CIT-EK-07A	•		•	•	•	•	•	•			10
CIT-EK-07B	•		•	•	•	•	•		•		10
CIT-EK-07C	•		•	•	•	•	•			•	10
CIT-EK-08A	•		•	•	•	•	•	•		•	10
CIT-EK-08B	•		•	•	•	•	•		•	•	10

Ref. Code	Sponge with Concentrated Enzyme	Sponge without Concentrated Enzyme	Eco Bowl 700ml	Transport Pad	Biohazard Label	Gauze Pads 10pcs	Pre-cleaning Valve Kit	4-piece Valve Kit	Lubricant Jelly, 42g	Lubricant Jelly, 20g	Kits /Box
CIT-EK-05	•		•	•			•	•			10
CIT-EK-05A		•	•	•			•	•			10