



**Quality First
Develop Harmoniously**

EN ISO 13485:2016

CE

FDA



CITEC™ ENDOSCOPY PROCEDURE KIT

CITEC™ Endoscopy Procedure Kit



Endoscopy units, as well as the health authorities in the world, for many years have been focusing their Infection control efforts on finding ways to minimize cross-contamination risks and improve patient care. As hospital acquired infections (HAIs) continue to pose a threat to patients and staff. **CITEC™ Endoscopy Procedure Kit** offer a solution for minimizing cross-contamination, improving patient care as well as improving the procedural efficiency.

CITEC™ Endoscopy Procedure Kit includes 9 disposable components, provides essential supplies needed for the procedure set up, bedside pre-cleaning, transport, and reprocessing of flexible endoscope.

Kit includes



Sponge with Concentrated Enzyme
Ergonomically round shape, highly match with the insertion part of endoscope. With multi enzyme detergents, to remove harmful bio-burden on the endoscopes.



Eco Bowl
Made with degradable material.
Maximum 700ml for detergent solution.



Transport Pad
Strong leak-proof liner to prevent leakage of fluids. Cinching cords converts pad into a transport container.




Biohazard Label
Stick the Biohazard Label on the transport pad, to warn the used endoscope.



Gauze Pads
Contain 10pcs, Size: 10*10cm.



Lubricant Jelly
20g or 42g, two kinds for selection.



Brush Kit
Double Heads Cleaning Brush: ϕ 5mm brush head*2, working length 2300mm, sheath ϕ 1.8mm.
Double Heads Valve Brush: ϕ 5mm brush head*1 and ϕ 10mm brush head *1.



Endoscope Valves System
For Olympus endoscopes series 140/160/180/190/240/260/290. Air Water Valve*1, Suction Valve*1, Biopsy Valve*1 and Water Jet Connector*1.



Pre-cleaning Valve Kit
For Olympus endoscopes series 140/160/180/190/240/260/290. AW Channel Cleaning Adapter*1, Suction Cleaning Adapter*1.

Specification:

| Ref. Code | Sponge with Concentrated Enzyme | Eco Bowl 700ml | Transport Pad | Biohazard Label | Gauze Pads 10pcs | Lubricant Jelly, 42g | Lubricant Jelly, 20g | Brush Kit | Endoscope Valves System | Pre-cleaning Valve Kit | Kits /Box |
|------------|---------------------------------|----------------|---------------|-----------------|------------------|----------------------|----------------------|-----------|-------------------------|------------------------|-----------|
| CIT-EK-02 | ● | ● | | | | | | | | | 5 |
| CIT-EK-04 | ● | ● | ● | ● | | | | | | | 5 |
| CIT-EK-06 | ● | ● | ● | ● | ● | | | ● | | | 5 |
| CIT-EK-07A | ● | ● | ● | ● | ● | ● | | ● | | | 5 |
| CIT-EK-07B | ● | ● | ● | ● | ● | | ● | ● | | | 5 |
| CIT-EK-08A | ● | ● | ● | ● | ● | ● | | ● | ● | | 5 |
| CIT-EK-08B | ● | ● | ● | ● | ● | | ● | ● | ● | | 5 |
| CIT-EK-09A | ● | ● | ● | ● | ● | ● | | ● | ● | ● | 5 |
| CIT-EK-09B | ● | ● | ● | ● | ● | | ● | ● | ● | ● | 5 |

