


Esophagus Anastomosis Stent



Application: To cover esophageal fistulas or in case of spontaneous or iatrogenic ruptures.

Performance Characteristics:

- Self expanding with high radial force
- Excellent positional stability
- Covered with PU/silicone, with good resistant and elastic performance
- Atraumatic ends
- Extraction threads on both ends for easy removal and repositioning
- 6 Radiopaque gold markers for better visibility



Specifications:

Product Code	Dia. A (mm)	Total Length B (mm)	Product Illustration	Introducer
NTJB4-AB-HC	28	80/100/120		18Fr*650mm Proximal Release
NTJB4-AB-FC				

Remark: HC means half Covered, FC means full Covered.