



Full covered Y Shape Tracheal/Bronchial Stent



Application: Y Shape Tracheal/Bronchial stent is used for tracheal and bronchial stenosis caused by various malignant lesions.

Performance Characteristics:

- Y-design fitting anatomic structure
- Self-expanding with high radial force
- Excellent positional stability
- Resistant and elastic covering
- Fully covered stents available
- Release with high radiopacity
- Golden radiopaque markers for clearer visibility



Specifications:

Product Code	Tracheal/ Bronchial Dia (mm)	Tracheal/ Bronchial right/ Bronchial left Length(mm)	Covering	Product Illustration	Introducer
NTJA2-LT1-16040	16/12	40/20/30	Half covered		24Fr*600mm
NTJA2-LT2-16050	16/12	50/20/30			
NTJA2-LT3-16060	16/12	60/20/30			
NTJA2-LT4-18040	18/14	40/20/30			
NTJA2-LT5-18050	18/14	50/20/30			
NTJA2-LT6-18060	18/14	60/20/30			
NTJA2-LT7-20040	20/14	40/20/30			
NTJA2-LT8-20050	20/14	50/20/30			
NTJA2-LT9-20060	20/14	60/20/30			
NTJA2-LT10-16040	16/12	40/15/30	Fully covered		24Fr*600mm
NTJA2-LT11-16050	16/12	50/15/30			
NTJA2-LT12-16060	16/12	60/15/30			
NTJA2-LT13-18040	18/14	40/15/30			
NTJA2-LT14-18050	18/14	50/15/30			
NTJA2-LT15-18060	18/14	60/15/30			
NTJA2-LT16-20040	20/14	40/15/30			
NTJA2-LT17-20050	20/14	50/15/30			
NTJA2-LT18-20060	20/14	60/15/30			